

TAPERED SHAFT™ FLOW CONTROLLER

WHEN TOTAL PUNCTAL OCCLUSION IS TOO MUCH

Eagle Vision's Tapered Shaft Flow Controller is designed for regulated flow of lacrimal fluid and maximum retention. The flesh-like comfortable flow controller is easy to insert and is reversible.



CE Marked

- Regulated flow of lacrimal fluid
- Maximum retention
- Easy insertion
- Flesh-like comfort
- Reversibility

Supplied	Description	REF
1 each single sterile silicone flow controller per pack, supplied pre-loaded, on a disposable inserter/dilator, color coded according to size	0.5mm diameter 0.6mm diameter 0.7mm diameter 0.8mm diameter	3065 3066 3067 3068

U.S. Patent Nos. 4,915,684 & 5,334,137 & 5,741,292



8500 Wolf Lake Dr., Suite 110 • P.O. Box 34877 • Memphis, TN 38133 USA
1-800-222-7584 • (901) 380-7000 • Fax (901) 380-7001
www.eaglevis.com • info@eaglevis.com