

# Coroneo™ Punctal Gauge\*



The Eagle Vision Coroneo Punctal Gauge is a complete gauging/dilating system in one instrument. It can be used to both gauge and dilate in increments from 0.5mm up to 1.0mm diameters. One tip is demarcated for 0.5, 0.6, and 0.7; flip it over for 0.8, 0.9, 1.0 on the other tip. So complete it takes the place of multiple instrument gauge sets.



CE Marked

- Single instrument gauge
- No need for multiple instruments
- Graduated demarcated tips on both ends
- Can be used with all punctum plugs
- Gauges and dilates from 0.5mm to 1.0mm diameters

Supplied	Description	REF
1 Gauge per package	Coroneo Punctal Gauge	2040

\*Developed in conjunction with Minas Coroneo, M.D.,  
New South Wales, Australia

Patents Pending



8500 Wolf Lake Dr., Suite 110 • P. O. Box 34877 • Memphis, TN 38133 USA  
1-800-222-7584 • (901) 380-7000 • Fax (901) 380-7001  
www.eaglevis.com • info@eaglevis.com