

# EaglePlug® TearFlow™

**Up to a  
150% larger lumen  
than our regular  
Flow Controller\***

Eagle Vision's patented flow controller design is now available with a larger lumen. The larger lumen is designed to allow additional flow in cases where total punctal occlusion is not desired.



CE Marked

- Designed for increased tear flow when compared to our original Tapered-Shaft™ Flow Controller
- Patented low profile rim for comfort
- Patented Eagle Vision tapered shaft for retention
- Sizes with larger lumen for modulated tear flow
- Easily reversible

PLUG DIAMETER	LUMEN DIAMETER	REF
0.5mm	0.2mm	3005
0.6mm	0.3mm	3006
0.7mm	0.3mm	3007
0.8mm	0.3mm	3008
0.9mm	0.4mm	3009
<b>NEW!</b> 1.0mm*	0.5mm	3010

\*When compared to our traditional Flow Controller with a 0.2mm lumen.

Supplied: 1 ea. sterile silicone TearFlow per single pack, two per box, preloaded on a disposable inserter.

U.S. Patent Nos. 5,334,137 & 6,027,470 and other Patents Pending.



8500 Wolf Lake Dr., Suite 110 • P.O. Box 34877 • Memphis, TN 38133 USA  
1-800-222-7584 • (901) 380-7000 • Fax (901) 380-7001  
www.eaglevis.com • info@eaglevis.com